

DIGITAL IFT-UST COMBO

an unique System for Physiotherapy Applications



FEATURES

- Usage of two modalities with individual timer is possible
- Feather-touch keypad
- 175 therapeutic-library programs
- Graphical LCD display (240 x 128)
- Individual ON-OFF control for all modalities
- SMPS powered
- Equipped with "ARM" microcontroller
- Cable check facility
- 30 configurable Favorites memory location

SPECIFICATIONS

IFT

- Mode : 4P, 4PV, 2P and Russian Current
- Frequency : 2, 3, 4, 5, 6, 7, 8, 9 and 10 KHz
- Base : 1 to 150Hz, 1 Hz Step
- Spectrum : 1 to 150Hz, 1 Hz Step, LIN, NL and TRP shape
- Vector type : 2PL, 4PL, 4PT & 4PNL
- Output : Isolated 0 – 98mA

UST

- Frequency : 1 MHz
- Mode : Continuous and pulse modes 1:1, 1:2, ... to 1:16
- Power : 0 to 3.1 W/cm²
- Timer : 1 – 30 minutes adjustable
- Input Power : SMPS @ 150V – 275V AC, 50 Hz Operated
- Dimension : 310 x 220 x 90 mm
- Weight : 2.9 Kg



No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics