

SEMI DIGITAL 4 CHANNEL TENS

an unique System for Physiotherapy Applications



FEATURES

- 3 modes of operation
- SMPS powered
- Stand alone Start/ Stop switch for output control
- Treatment cycle and other parameters are configurable
- Equipped with "ARM" microcontroller
- User friendly 7 segment display and Feather-touch Keys as user interface

SPECIFICATIONS

- Mode : Continuous, Burst and Sweep
- Channels : Four
- Frequency : 1 to 250Hz @ 400uS width
- Sweep : 9 programs, SP0, SP1, SP2, SP3, SP4, SP5, SP6, SP7, SP8 and SP9
- Burst : 3 Sec Hold & 2 Sec Rest
- Output : 0-99mA adjustable and isolated
- Timer : 1 – 30 minutes adjustable
- Input Power : SMPS @ 100V – 275V AC, 50 Hz Operated
- Dimension : 270 x 250 x 90 mm
- Weight : 2 kg



No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics