

SEMI DIGITAL DUAL UST

an unique System for Physiotherapy Applications



FEATURES

- Ruggedized probe holder
- SMPS powered
- Stand alone Start/ Stop switch for output control
- Treatment cycle and other parameters are configurable
- Equipped with "ARM" microcontroller
- User friendly 7 segment display, Bar Graph LED and Feather-touch Keys as user interface
- Isolated probe with elevated safety
- Single probe for both 1 and 3 MHz output

SPECIFICATIONS

- Frequency : 1 and 3 MHz
- Mode : Continuous & Pulse
- Pulse Ratio : 1:1,1:2, 1:4, 1:7,1:8 ,1:12 & 1:16
- Power : 0.1 to 3.1 watts/cm²

- Timer : 1 – 30 minutes adjustable
- Input Power : SMPS @ 100V – 275V AC, 50 Hz Operated
- Dimension : 270 x 250 x 90 mm
- Weight : 2.2 Kg

No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083

+91-44-4775 9102

+91-73586 20147 / 73586 20247

alphameditronics@gmail.com

www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics