

SEMI DIGITAL IFT

an unique System for Physiotherapy Applications



FEATURES

- Isolated dual channel Output
- SMPS powered
- Stand alone Start/ Stop switch for output control
- Treatment cycle and other parameters are configurable
- Equipped with "ARM" microcontroller
- User friendly 7 segment display, Bar Graph LED and Feather-touch Keys as user interface

SPECIFICATIONS

- Mode : 4P, 4PV-2PL, 4PV-LIN, 4PV-TRP, 4PV-NL, 2P and Russian Current
- Frequency : 2, 3, 4, 5, 6, 7, 8 and 9 KHz
- Base : 1 to 250Hz, 1 Hz Step
- Spectrum : 0 to 99Hz, 1 Hz Step, LIN shape
- Output : Isolated 0 – 97mA
- Timer : 1 – 30 minutes adjustable
- Input Power : SMPS @ 100V – 275V AC, 50 Hz Operated
- Dimension : 270 x 250 x 90 mm
- Weight : 2.3 kg

No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083

+91-44-4775 9102

+91-73586 20147 / 73586 20247

alphameditronics@gmail.com

www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics