

SEMI DIGITAL LASER

an unique System for Physiotherapy Applications



FEATURES

- Ruggedized probe holder
- SMPS powered
- Stand alone Start/ Stop switch for output control
- Treatment cycle and other parameters are configurable
- Equipped with "ARM" microcontroller
- User friendly 7 segment display, Bar Graph LED and Feather-touch Keys as user interface
- Single therapy probe for IR and RED

SPECIFICATIONS

- Mode : IR, RED and IR + RED
- IR wavelength : 780 nm
- RED wavelength : 660 nm
- Energy Level : 0.1 to 25 J/cm²
- Frequency : 1 to 9.999 KHz , 500 Hz step size

- Timer : 1 – 8 minutes adjustable
- Input Power : SMPS @ 100V – 275V AC, 50 Hz Operated
- Dimension : 270 x 250 x 90 mm
- Weight : 2 Kg



No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics