

TENS Plus

an unique System for Physiotherapy Applications



FEATURES

- Isolated and independent two channel output channels
- Configurable frequency and intensity level using
- 7 segment display and Feather-touch Keys
- Equipped with “ARM” microcontroller
- Stand alone Start/ Stop switch for output control
- Built-in 15 minutes timer
- 10 Sweep programs

SPECIFICATIONS

- Mode : Continuous, Burst and 10 sweep modes
- Frequency : 1 to 250 Hz @ 400 μ s ON time
- Sweep program : 10 programs, SP0,SP1,SP2,SP3,SP4,SP5,SP6,SP7,SP8 and Sp9
- Burst period : Hold @ 3 sec, Rest @ 2 Sec
- Timer : in built 15 minutes
- Output : 0 – 99 mA adjustable
- Power Supply : 12V Adaptor Operated
- Dimension : 90 x 133 x 35 mm
- Weight : 350 g

No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083

+91-44-4775 9102

+91-73586 20147 / 73586 20247

alphameditronics@gmail.com

www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics